

# PCS 7000 Oxidizing Disinfectant/Disinfectant Cleaner

## Description:

Disinfectant/Disinfectant Cleaner for use on hard non-porous environmental surfaces in health care facilities, institutions, schools and hospitality industries, where organic soils may be present. Use to kill or inactivate spores of *C. difficile* in settings where contamination by fecal matter is likely.

- Kills Bacteria (Bactericide) • Kills *Pseudomonas aeruginosa*, *Salmonella enterica*, and *Staphylococcus aureus*
- Kills *Clostridium difficile* (*C. difficile*) (spore form)

Guarantee 0.70% Sodium Hypochlorite when packed.  
DIN: 02314878

## Directions for use.

Pull up tab and remove large cover. Thread first wipe in centre of roll through slits in the cover. DO NOT PUSH FINGER THROUGH OPENING Replace cover. Pull wipe up and slightly to the side. Always snap lid securely shut between uses to prevent moisture loss. To disinfect walls, tables, counter tops. Use a fresh sheet to wipe pre-cleaned surfaces. Wipe surface and keep wet for 5 minutes. Rinse with potable water all surfaces and/or objects that come in contact with food or objects that come in contact with children at the mouthing stage of development. Avoid contamination of food.

## Product Codes:

- #6063-6 • 6 x 946 mL w/ flip top
- #6068-2 • 4 x 2.5 L
- #6030-2 • 4 x 3.78 L closed loop for dispenser use.
- #6030-4 • 2 x 3.78 L

